

# Cell-based NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling

## NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling in HeLa cells: Phosphorylation and degradation

### Introduction

Inactive NF- $\kappa$ B is located in the cytosol in a complex with the inhibitory protein I- $\kappa$ B $\alpha$ . A variety of extracellular stress-related signals, including inflammatory cytokines, free radicals, UV irradiation, and bacterial or viral antigens can initiate phosphorylation and activation of I- $\kappa$ B kinase (IKK). Activated IKK phosphorylates I- $\kappa$ B $\alpha$ , which subsequently dissociates from NF- $\kappa$ B. Phosphorylated I- $\kappa$ B $\alpha$  is ubiquitinated and degraded by the proteasome. Phosphorylated NF- $\kappa$ B is then translocated into the nucleus where it acts as a transcriptional activator.

The NF- $\kappa$ B pathway plays an important role in regulating immune responses to infection, and deregulation of NF- $\kappa$ B has been linked to several diseases, including cancer, and autoimmune diseases.

### Materials used


HeLa cells are available from ATCC. Cell culture media and additives are available from Invitrogen. TNF $\alpha$  and TPCA-1 are available from Sigma. AlphaScreen SureFire phospho NF- $\kappa$ B, total NF- $\kappa$ B, phospho I- $\kappa$ B $\alpha$ , and total I- $\kappa$ B $\alpha$  assay kits are available from PerkinElmer. Western blotting detection reagents and secondary antibody-HRP conjugates are available from Perkin Elmer. Tissue culture plasticware is from Nunc, and AlphaScreen SureFire assay plates are available from PerkinElmer.

### Results

#### Stimulation time point selection

HeLa cells were stimulated with 10 ng/mL TNF $\alpha$  for various lengths of time at room temperature. The cells were lysed, and the lysates were analysed for phosphorylated and total levels of both NF- $\kappa$ B (Figure 1) and I- $\kappa$ B $\alpha$  (Figure 2), in order to determine the optimal time point for measuring phosphorylation of both NF- $\kappa$ B and I- $\kappa$ B $\alpha$ .

Under the conditions tested, maximal phosphorylation of NF- $\kappa$ B was observed after 5 minutes of stimulation. Elevated phosphorylation levels were observed until 10 minutes after which levels decreased. Total protein levels of NF- $\kappa$ B remained unchanged for the duration of the time course stimulation.




Similarly, I- $\kappa$ B $\alpha$  phosphorylation stimulated by TNF $\alpha$  was also transient. Maximal levels of phosphorylation were observed at 5 minutes, and returned to basal levels after 15 minutes. Furthermore, marked differences in the level of total I- $\kappa$ B $\alpha$  were observed. Levels of I- $\kappa$ B $\alpha$  rapidly decreased after 5 minutes of TNF $\alpha$  stimulation, and then gradually started to increase again after 30 min. I- $\kappa$ B $\alpha$  phosphorylation and degradation were confirmed by western blotting (Figure 2). To test if

# Cell-based NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling

## NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling in HeLa cells: Phosphorylation and degradation


the decrease in total I- $\kappa$ B $\alpha$  observed was a result of ubiquitination blocking the antibody recognition, a different antibody - raised to a distal epitope - was also tested in a Western blot. A similar depletion of total I- $\kappa$ B $\alpha$  result was observed (data not shown).



Based on these results, 5 minute stimulation with TNF $\alpha$  at room temperature was selected as a useful time point for further analysis of NF- $\kappa$ B and I- $\kappa$ B $\alpha$  phosphorylation.

### Dose-dependent stimulation with TNF $\alpha$

HeLa cells were stimulated for 5 minutes with various concentrations of TNF $\alpha$ , then lysed. The lysates were analysed for phosphorylated and total levels of both NF- $\kappa$ B and I- $\kappa$ B $\alpha$ . TNF $\alpha$ -induced phosphorylation of both proteins in a dose dependent manner with EC<sub>50</sub> values of 0.74 ng/mL and 2.4ng/mL for NF- $\kappa$ B and I- $\kappa$ B $\alpha$ , respectively (Figure 3).



# Cell-based NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling

## NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling in HeLa cells: Phosphorylation and degradation

Total levels of NF- $\kappa$ B protein remained unchanged as expected, at all concentrations of TNF $\alpha$  treatment. However, total levels of I- $\kappa$ B $\alpha$  protein were slightly reduced at higher concentrations of TNF $\alpha$ , suggesting that at higher concentrations of TNF $\alpha$ , I- $\kappa$ B $\alpha$  protein degradation may occur at a timepoint as early as 5 minutes after stimulation (Figure 4).




Figure 4

### Inhibition with TPCA-1

The effect of the I- $\kappa$ B kinase 2 (IKK $\beta$ ) selective inhibitor, TPCA-1, on NF- $\kappa$ B and I- $\kappa$ B $\alpha$  phosphorylation was examined. HeLa cells were treated with various concentrations of TPCA-1 for 60 minutes, and then stimulated with 10 ng/mL TNF $\alpha$  for 5 minutes. The cells were lysed, and the lysates were analysed for modulation of NF- $\kappa$ B and I- $\kappa$ B $\alpha$  (Figure 5). TPCA-1 was found to inhibit phosphorylation of both targets in a dose dependent manner with IC<sub>50</sub> values of 5.8  $\mu$ M and 0.67  $\mu$ M for NF- $\kappa$ B and I- $\kappa$ B $\alpha$ , respectively.

Again, total levels of NF- $\kappa$ B protein were unchanged. Total levels of I- $\kappa$ B $\alpha$  protein were slightly reduced at lower concentrations of inhibitor (data not shown).




Figure 5


### I- $\kappa$ B $\alpha$ protein degradation

In the previous experiments, 5 min stimulation with TNF $\alpha$  was used, as this time point resulted in maximal phosphorylation of both proteins. As the NF- $\kappa$ B/I- $\kappa$ B $\alpha$  pathway is regulated by degradation of I- $\kappa$ B $\alpha$  as well as phosphorylation events, the dose response experiments were repeated, stimulating for 15 minutes with TNF $\alpha$  rather than 5 minutes, to ensure maximal loss of I- $\kappa$ B $\alpha$  signal. TNF $\alpha$  treatment for this period of time induced dose-dependant I- $\kappa$ B $\alpha$  degradation (EC<sub>50</sub> 0.34 ng/mL), doses comparable to the doses observed for NF- $\kappa$ B phosphorylation (Figure 6). TNF $\alpha$  induced I- $\kappa$ B $\alpha$  protein degradation could be inhibited by pre-treatment with TPCA-1, with an IC<sub>50</sub> of 1.3  $\mu$ M, a value similar to that observed from direct

# Cell-based NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling

## NF- $\kappa$ B/I- $\kappa$ B $\alpha$ signaling in HeLa cells: Phosphorylation and degradation

inhibition of phosphorylation of I- $\kappa$ B $\alpha$  measured at 5 min (Figure 5).



### Conclusions

Modulation of the NF- $\kappa$ B/I- $\kappa$ B $\alpha$  pathway is readily detectable in HeLa cells in response to TNF $\alpha$  treatment. Using AlphaScreen SureFire assays, phosphorylated NF- $\kappa$ B and I- $\kappa$ B $\alpha$  can be detected. For I- $\kappa$ B $\alpha$ , both phosphorylation and also total protein levels could be used as a marker for pathway regulation. In contrast, NF- $\kappa$ B levels were not modulated by TNF $\alpha$ , and phosphorylation was the most useful marker.

### Ordering Information

All cell lines are available from ATCC ([www.atcc.org](http://www.atcc.org)).

AlphaScreen<sup>®</sup> SureFire<sup>®</sup> cellular assay kits, available in 500pt, 10,000pt and 50,000pt pack sizes, are from PerkinElmer.

Catalog#

Note: These kits must be used in conjunction with the AlphaScreen protein A