


# TF-1 cells

## Handling suggestions for TF-1 cells

The TF-1 cell line was established from a heparinized bone marrow aspiration sample from a 35 year old male with severe pancytopenia. The cells are completely dependent on IL-3 or GM-CSF for culturing, but in the presence of these factors grows rapidly. TF-1 cells respond to a variety of cytokines, and is a useful tool for studying cytokine signaling mechanisms because they carry a wide range of endogenous receptors. Phosphorylation of endogenous STAT1, STAT3 and STAT5 is readily detectable in these cells, and provides a convenient marker for cytokine response.



Timecourse of stimulation of STAT phosphorylation in TF-1 cells. STAT 1 phosphorylation is stimulated by IFN $\gamma$ , STAT3 phosphorylation is stimulated by IL-6 and STAT5 phosphorylation is stimulated by GM-CSF. Phosphorylation detected by AlphaScreen<sup>®</sup> SureFire<sup>®</sup> cellular assay kits.

**Morphology:** lymphoblast

**Source:** blood

**Growth:** suspension

**Organism:** Human

**Sources:** ATCC: CRL-2003

ECACC Cat# 93022307

**Suggested media:**

RPMI (Gibco Cat#11885) supplemented with 10% FBS (Gibco) and 2 ng/mL GM-CSF (R&D Systems, Cat# 215-GM)

1% Sodium pyruvate (Gibco Cat#11360)

1% Pen Strep Glutamine (Gibco Cat#10378)

**Culturing suggestions:**

Maintain cells at densities in the range of  $5 \times 10^4$  –  $5 \times 10^5$  cells/mL.

Add fresh media as required.

**Detectable signaling pathways:**

STAT

**Known Receptors:**

IL-3, IL-6, GM-CSF, IFN $\gamma$